

NASA/Goddard
10 - 11 Feb. 1999

CsI Crystal Testing: Wrappings and End Treatments

J. Eric Grove
Naval Research Lab
11 February 1999




Position Response of CsI

NASA/Goddard
10 - 11 Feb. 1999

Map of scintillation in 36 cm CsI crystal.

- ❑ Scanned ^{137}Cs source (662 keV).
- ❑ Crystal viewed full-face by PMT (connected with optical grease). Far face was blackened.
- ❑ Side wrap was Tetrtek and aluminized mylar.
- ❑ Scintillation light yield drops by \sim half over length of crystal.
- ❑ Solid line is 42-cm exponential attenuation length.
- ❑ “Hotspot” at \sim 13 cm is real. \sim 2-3% magnitude similar to BaBar hotspots.




Long-term Pressure Tests

NASA/Goddard
10 - 11 Feb. 1999

36-cm crystal, Tetrtek and mylar wrap,
held under pressure and scanned.

- ❑ >50 days of 10-g load on all surfaces.
 - All curves normalized to first measurement.
- ❑ Light yield decreases under pressure.
- ❑ Light yield stabilizes after ~10 days at ~5 - 15% loss.
- ❑ Pressure increased to 23 g.
- ❑ Light yield rapidly stabilizes at an additional ~3% loss.




Naval Research Lab
Washington DC



Crystal Test Procedure

NASA/Goddard
10 - 11 Feb. 1999

- ❑ Crystals are numbered and inspected as delivered.
 - Factory wraps are teflon-only or Tyvek and aluminum foil.
 - Additional wrap of aluminized mylar is added if necessary.
- ❑ Various end treatments may be applied as required for test.
- ❑ Crystal is mounted above 2" PMT.
 - Air gap between crystal and PMT.
- ❑ ^{22}Na or ^{137}Cs source is scanned along crystal with motor drive.
- ❑ **Good spectroscopy is achieved.**



Naval Research Lab
Washington DC





Ukrainian Crystal Tests


NASA/Goddard
10 - 11 Feb. 1999

- ❑ Crystals received from Ukraine
 - Four 310 x 30 x 23 mm
 - Two 400 x 30 x 23 mm

- ❑ Test data accumulated at Ukraine.
 - Crystal viewed by PMT with air gap.
 - Surfaces polished and roughened.
 - Wrap is Tyvek and aluminum foil.
 - Far end is Tyvek.
 - ^{22}Na source scanned along length.

γ -Na-22, $E_\gamma=1275\text{keV}$

167P3.31.17




White Mask

NASA/Goddard
10 - 11 Feb. 1999

32-cm crystal scanned with ^{22}Na .

- ❑ All surfaces polished. Tetratex wrap.
- ❑ Viewed by PMT with air gap.
- ❑ Near face masked with Tyvek.
- ❑ Two masks, different apertures:
 - Size and location of large PIN.
 - Size and location of small PIN.
- ❑ Light tapering is independent of aperture size.
- ❑ Attenuation length (beyond 10 cm)

$$\lambda = 110 \text{ cm.}$$



Naval Research Lab
Washington DC




Black Mask

NASA/Goddard
10 - 11 Feb. 1999

32-cm crystal scanned with ^{22}Na .

- ❑ All surfaces polished. Tetratek wrap.
- ❑ Viewed by PMT with air gap.
- ❑ Near face masked with black paper.
- ❑ Two masks, different apertures:
 - Size and location of large PIN.
 - **Size and location of small PIN.**
- ❑ **Light tapering is independent of aperture size.**
- ❑ Attenuation length (all crystal)
 $\lambda = 75$ cm.




Naval Research Lab
Washington DC



Black or White Ends?


NASA/Goddard
10 - 11 Feb. 1999

- ❑ How does end treatment affect light yield and attenuation?
 - 32-cm Crismatec crystal mapped with ^{22}Na source.
 - PMT readout with black or white aperture mask ($\sim 1 \text{ cm}^2$ open).
- ❑ Black mask reduces light S to $\sim 2/3$ of white mask.
- ❑ Black mask shortens attenuation length.
 - $\lambda = 75 \text{ cm}$ for black
 - $\lambda = 110 \text{ cm}$ for white
- ❑ Position resolution scales as λ / \sqrt{S}
 - **Black mask gives**
 - 1/3 less light, but
 - 20% better position resolution.



Naval Research Lab
Washington DC






Two-ended Measurements


NASA/Goddard
10 - 11 Feb. 1999

- ❑ View crystal from both ends through black mask with aperture of size and location of large PIN.
 - “Red” = Hamamatsu R669 PMT
 - “Blue” = Hamamatsu R5900 PMT
- ❑ Both ends show same attenuation length
 $\lambda = 70$ cm
- ❑ ooops. Vibration knocked blue PMT out of fixture during run.



Naval Research Lab
Washington DC






Light Asymmetry

NASA/Goddard
10 - 11 Feb. 1999


- ❑ We use light asymmetry to measure position of shower (or photopeak) along crystal.
- ❑ Light asymmetry measure:
 - If attenuation is linear, use
 - $x = 2 \lambda (\text{Right} - \text{Left}) / (\text{Right} + \text{Left})$
 - If attenuation is exponential, use
 - $x = \lambda \log(\text{Right} / \text{Left})$
- ❑ Note that this crystal has twice the slope of the 32-cm crystal used at SLAC and CERN.
 - Will give factor-of-two better position resolution.
- ❑ Asymmetry can be readily mapped with lab sources and PMTs.



Naval Research Lab
Washington DC



- ❑ We mapped the light yield across the face of the crystal.
 - Varied size and location of aperture across face of crystal.
 - Aperture sized for big PIN at different locations on end face.
- ❑ The light yield varies across the end face of the crystal.
 - Center and corners are bright.
 - Sides are dim.
 - Note total effect is only 20-30%.



CERN Beam Test '98

NASA/Goddard
10 - 11 Feb. 1999

